

Equipment for A.I. stations of bulls

Collection of semen

Artificial vagina – rubber roller, valve
Inner of artificial vagina with notch carving
Plastic collector of semen



Optidyl – concentrated medium ready to use for bull semen

Tris formula well known all over the world

Content: tris diluent, ionised egg yolk, glycerol, antibiotics
Antibiotic formula according to EC requirements: penicillin,
streptomycin, spectinomycin, lincomycin

preparation :

mix the content of one flask of 500 ml with 750 ml of bidistilled
water

storage at +5°C




Cryo Box for manipulation with A.I. doses in liquid nitrogen




Guiding rod for bull's



Straws 0.25 and 0.5 ml



Straws 0,25 ml / straws 0,5 ml



AM160 / AN160 white

AM163 / AN163 clear

AM261 / AN261 blue roi

AM262 / AN262 clear blue

AM264 / AN264 turquoise blue

AM265 / AN265 pale blue

AM360 / AN360 fluo yellow

AM361 / AN361 sun flower yellow

AM362 / AN362 pale yellow

AM363 / AN363 clear yellow

AM460 / AN460 dark pink

















AM461 / AN461 salmon

AM463 / AN463 red

AM465 / AN465 dark red

AM466 / AN466 clear red

AM468 / AN468 fluo pink

	AM560 / AN560 pistachio green
	AM561 / AN561 green
	AM563 / AN563 fluo green
	AM566 / AN566 clear pale green
	AM567 / AN567 dark green
	AM660 / AN660 dark purple
	AM661 / AN661 pale purple
	AM662 / AN662 clear purple
	AM761 / AN761 fluo orange
	AM763 / AN763 clear orange
	AM765 / AN765 dark orange
	AM960 / AN960 putty
	AM961 / AN961 maroon
	AM963 / AN963 grey
	AM964 / AN964 silver grey
	AM965 / AN965 black

Packing

2000 pcs Bag

Automatic filling sealing printing machine

BOVINE ARTIFICIAL INSEMINATION

CRYO-VET
Artificial insemination equipment
PIG'INS

MX4
Fully automatic machine
for 0.25 and 0.50 ml
semen straws

MX4
Filling, sealing
and printing machine
for the straws

Patent pending

MX4
From 1 to 4 straws are filled
and sealed per cycle

Main advantage:
0% SEMEN LOST IN THE PLUG

MX4
15 000 straws per hour.

